

**skin condition**

Balanced to mature or prematurely-aging skin.

**description**

Pro Microneedling is an advanced treatment that delivers innovative techniques and powerful, professional-grade actives to brighten the complexion, refine the look of pores, reduce the appearance of dark spots, and diminish the appearance of lines and wrinkles.

**benefits**

- Brightens complexion
- Reduces the appearance of dark spots
- Diminishes the appearance of lines and wrinkles
- Refines the look of pores

**featured professional products**

- Pro Power Peel
- Pro Restore
- Exfoliant Accelerator 35



# pro microneedling

**how it works****1. Preparation**

Skin resurfacing with customized Pro Power Peel

**2. Treatment**

Microneedling with professional serums

**3. Recovery**

In-treatment room recovery along with post-treatment care and home care

**treatment protocol****double cleanse**

1. PreCleanse
2. Perform Face Mapping
3. Skin specific cleanser

**resurfacing**

1. Apply One Step Prep.
2. Apply specific Pro Power Peel (1-2 layers).
3. Apply Neutralizing Solution after Pro Power Peel and remove as directed.
  - a. Power Down: Use Exfoliant Accelerator 35 or omit resurfacing.

- Dull Skin: BioLumin-C Pro Serum or Pro Restore at 0.2-0.6mm
- Pore Refining: Niacinamide IonActive at 0.4-0.8mm or Pro Restore 1.0-1.2mm
- Hyperpigmentation: PowerBright IonActive Serum or Pro Restore at 0.2-0.8mm
- Lines and Wrinkles: Pro Restore at 0.4-1.2mm

**dermal layering**

1. Apply Pro Restore.

**deep treatment**

1. Apply Pro Restore
2. Perform the Dermalogica Microneedling Technique
  - Wipe the skin the damp esthetic wipes after microneedling
  - If including neck and chest, add additional 10 minutes per area.